


Example Evacuation Set-up



Fast and Easy Full Flow Evacuation

The Megaflow™ Basics Kits help technicians remove restrictive valve cores and perform AC/R system evacuations up to 10x faster. Maximize the full flow potential of any vacuum pump!


VCRT
VALVE CORE REMOVAL TOOL™

MegaFlow
VACUUM-RATED HOSES™



With both VCRTs, you can remove the cores from **1/4 in. and 5/16 in. Schrader Valves** with ease.

Valve cores are highly restrictive, **blocking up to 90% of all flow**. Not removing valve cores limits the maximum flow of any vacuum pump to **≈0.2 cfm**.




Each MegaFlow vacuum hose features a **1/2 in. internal diameter** for maximum flow and 1/4 in. - 3/8 in. fittings.

Multiple inputs and large-diameter hoses allow a vacuum pump to operate at **full capacity**. Hose diameter and quality will dramatically impact evacuation times!

Evacuating a 10-ton system to 500 microns

6-ft, 1/4" hose through valve core	40 min.
6-ft, 1/4" hose & core removed	19 min.
6-ft, 1/2" hose & core removed	3 min.

With these high conductance speed AC/R evacuation essentials, your vacuum pump can finally reach its true potential. The entire setup is vacuum-rated to 20 microns and designed for high-pressure refrigerants, including R-410A and R-454B. Full flow and faster evacuation guaranteed!



MegaFlow

BASICS™



Faster Evacuation Made Simple

Kit Includes:

MegaFlow
1/4 IN. BASICS™



Part# MGABAS

- (1) 6 ft. 1/2 in. Evacuation Hose
- (1) 1/4 in. Valve Core Removal Tool

MegaFlow
5/16 IN. BASICS™



Part# MGABRO

- (1) 6 ft. 1/2 in. Evacuation Hose
- (1) 5/16 in. Valve Core Removal Tool



SCAN TO LEARN MORE



SCAN TO LEARN MORE